

BIDAC
801-585-0968

Cell Imaging
801-587-7964

DNA Peptide
801-581-4051

DNA Sequencing
801-585-2976

Drug Discovery
801-585-1527

Electron
Microscopy
801-585-1242

Flow Cytometry
801-581-7004

Genomics
801-585-2977

Iron & Hematology
801-581-6870

Machine Shop
801-581-3180

Mass Spectrometry
801-581-5018

Material Sciences-
Engineering
801-581-5303

Metabolic
Phenotyping
801-587-7696

Metabolomics
801-587-7779

Mutation
Generation
801-585-0662

Nuclear
Engineering
801-587-9696

Nuclear Magnetic
Resonance
801-581-8418

Protein Interactions
801-585-5021

Recombinant
Protein
801-581-5311

Scalable Analytics
& Informatics
801-581-8067

Small Animal
Imaging
801-587-8357

Small Animal
Ultrasound
801-587-7849

Zebrafish
801-585-6574

Flow Cytometry Shared Resource Laboratory

Wintrobe Bldg. Rom 603
26 N 1900E
Salt Lake City, UT 84132
801-585-7382

ASSAY REPORT: FRESH VS FROZEN UMBILICAL CORD TISSUE AND CELL SAMPLE

Date: January 18-2017
Report Created by: JEM

Cell Viability and Identity Protocol

- Pre-warmed water bath to 37°
- Add sample ID to 15ml conical tubes. One tube for each frozen vial.
- Add 10ml wash buffer to each 15ml conical
- Put frozen vial into holes in floater and put into water bath
- As soon as vial is 50% thawed, take out of water bath and dump into 15ml conical tube.
- With pipette, transfer the remaining liquid from vial into conical tube.
- Spin conical tubes 250g X 3min.
- Aspirate all liquid from tubes.
- Add 1ml of wash buffer and gently re-suspend with pipette.
- Transfer 300ul to pre-labelled tubes
- Add individual antibodies or premade cocktail to appropriate tubes.
- Tube1=unstained
- Tube 2= Isotype control
- Tube 3= Sample
- Incubate 30min 4°
- Add 4ml wash buffer
- Spin 250g X 3min
- Re-suspend in 500ul wash buffer
- Add 10ul of Accudrop Counting beads to Tube 3
- Run on flow cytometer immediately.
- Add 5ul DAPI 5 min prior to acquisition

Cat#	Vendor	Antibody	Fluor	ul/test
11-0909-41	ebioscience	CD 90	FITC	2.5ul
11-4714-41	ebioscience	IgG	FITC	2.5ul
12-1057-41	ebioscience	CD105	PE	2.5ul
12-4714-41	ebioscience	IgG	PE	5ul
25-9459-41	ebioscience	CD45	PECY7	2.5ul
25-4714-41	ebioscience	IgG	PECY7	2.5ul
17-0739-41	ebioscience	CD73	APC	2.5ul
17-4714-41	ebioscience	IgG	APC	0.625
		viability	DAPI	5ul

Reagents

Wash Buffer- 2.5%BSA in PBS (Core)

ASSAY RESULTS

Fresh Sample ID

Excluding RBC contamination, cellularity consisted of 9.7×10^5 viable cells per mL. Viability based on exclusion of DAPI equals 67.6%. Cellular phenotype consisted of:
 CD90= 59.2%
 CD105= 55.4%
 CD73= 86.7%
 CD45= 3.19%

Total Viable Mesenchymal Stem Cell Count: 1,200,000 per mL. (This number is an approximation based on the use of bead markers. There is to be expected normal cell loss based on the steps to process the sample before Flow Cytometry testing. The 1,200,000 per mL number is the cell count that remains after testing protocol preparations.)

Frozen Sample ID

Excluding RBC contamination, cellularity consisted of 7.8×10^5 viable cells per mL. Viability based on exclusion of DAPI equals 43.1%. Cellular phenotype consisted of:
 CD90= 28.0%
 CD105= 35.1%
 CD73= 84.4%
 CD45= 4.06%

Total Viable Mesenchymal Stem Cell Count: 800,000 per mL. (This number is an approximation based on the use of bead markers. There is to be expected normal cell loss based on the steps to process the sample before Flow Cytometry testing. The 800,000 per mL number is the cell count that remains after testing protocol preparations.)

James Marvin
 Director, Flow Cytometry SRL